



SKeen

Original user manual

Software

skeen - Software

Original user manual

Please read the original user manual carefully before using the product.



Copyright © 2024 by FotoFinder Systems GmbH

Softwareversion from 1.0.2
Date: 19.12.2024



FotoFinder Systems GmbH
Industriestraße 12
84364 Bad Birnbach
Germany

www.fotofinder.de
info@fotofinder.de
Tel.: +49 (0) 8563 – 97720-0
Fax: +49 (0) 8563 – 97720-10

Content

1	Safety.....	5
1.1	Intended use.....	5
1.2	User groups	6
1.3	Patient population	6
1.4	Indications and contraindications.....	6
1.5	Foreseeable misuse.....	7
1.6	Residual risks.....	7
1.6.1	IT-Security.....	7
2	FotoFinder skreen software	8
2.1	First login	8
2.2	Installation und updates.....	8
2.3	PIN (Personal Identification Number)	8
2.4	FotoFinder Hub login	9
2.5	General operating information.....	10
2.5.1.1	Sleep mode	10
2.5.1.2	Software navigation	10
2.6	Home screen.....	10
2.7	Menu Bar.....	11
2.8	Patients.....	12
2.9	Search and select existing patients	12
2.9.1	Create new patient	13
2.9.2	Remove current patient from selection	13
2.10	Creating images.....	14
2.10.1	Saving the localisation	16
2.10.2	Saving the localisation at a later stage	16
2.11	Sessions	17
2.11.1	Assigning a session to another patient afterwards	17
2.12	The AI Score	18
2.12.1	Requesting the AI-Score (AIMEE).....	18
2.12.2	AI Score result information.....	19
2.13	About FotoFinder	20
2.14	Settings.....	21
3	Appendix	22

1 Safety

1.1 Intended use

The FotoFinder software works in conjunction with the FotoFinder Hub online cloud. The software is designed for patient management, standardized documentation of microscopic and macroscopic images, and to assist in the initial assessment of skin conditions. The FotoFinder software enables digital documentation of intact human skin by healthcare professionals. The microscopic images are stored together with the relevant patient data, which makes it possible to visualize changes in lesions during subsequent follow-up examinations of the patient.

The FotoFinder software is used in combination with FotoFinder or DermLite imaging devices, which allow to capture microscopic images using a mobile device.

The following features are available:

- Acquisition and management of patient data
- Capturing and managing microscopic and macroscopic images
- Documentation of images in sessions
- Assigning a session to a patient
- Assigning a body site (localization) to an image
- Request AI score (Artificial Intelligence)

The FotoFinder software connects online with the Moleanalyzer pro algorithms to generate the AI score. This function requires a specific license and plan, which are managed in the FotoFinder Hub account. The data is synchronized, stored and managed via this cloud solution.

The FotoFinder software is intended for skin lesions. The software must not be used to make or confirm a clinical diagnosis of melanoma, any other skin disease or skin cancer.

The software may only be used by physicians or healthcare professionals trained in the clinical diagnosis of skin cancer or other skin diseases.

The software does not provide a diagnosis. The AI score is based on statistics. The diagnosis and therapy decision are the responsibility of the physician!

The software is intended for transient use. In combination with the imaging device, the product is in continuous use for less than 60 minutes during a diagnosis session.

1.2 User groups

The following target groups with necessary qualifications may work with the application:

Target group	Qualification	Permitted work in the application
Medical or healthcare professionals	Trained and instructed, as well as professionally qualified by a completed professional training in the medical field.	<ul style="list-style-type: none"> - Accessing and managing patient records and respective data - Capturing images of skin - Viewing skin images - Diagnosing skin diseases - Treatment and therapy of skin and hair diseases

1.3 Patient population

Patients with one of the following characterizations are intended to be examined with the software:

- General persons with skin lesions, moles
- Patients with multiple nevus syndrome
- People with high risk of skin cancer / family history of skin cancer

The intended patient population includes patients regardless of demographic factors (e.g. gender, age, profession), physical factors (e.g. weight, height, strength) or social, religious and cultural background.

1.4 Indications and contraindications

The FotoFinder software may only be used for the following indications:

- Primary skin lesions with a diameter of up to 8 mm.
- Lesions on intact skin (i.e., non-ulcerated and non-bleeding lesions).
- Lesions without scarring or fibrosis due to previous trauma.
- Lesions that are not in close proximity to psoriasis, eczema, acute sunburn, or similar skin conditions
- Lesions that are not located on specific body sites (i.e., not on the acra, genitalia, eyes, mucous membranes, orifices)

Note the following contraindications for images when requesting an AI Score:

The captured lesion

- must not be covered by hair.
- must not have any foreign bodies or markings within a radius of 30 mm (i.e. tattoos, paint markings of any kind).

The algorithm was trained with images of Fitzpatrick skin type I-III. Do not use the AI Score on patients with skin type IV or higher, as the performance of the algorithm was not assessed and therefore the accuracy of the algorithm cannot be claimed.

NOTE

Clinically unambiguous melanomas may not be evaluated for the purpose of making a treatment decision with the FotoFinder software.

1.5 Foreseeable misuse

The following points describe foreseeable misuse of the software:

- The physician incorrectly assumes that the software provides a diagnosis.
- The physician bases their diagnosis exclusively on results of software.
- The application for documentation is performed on non-intact skin, mucous membranes or in body orifices.
- The physician believes that the accuracy of the AI Score can be claimed and assumes that the score is indicative of the malignancy of the mole.
- The physician requests an AI Score for an image that does not meet the requirements, e.g., due to body hair, visible tattoo, or size of the lesion.

NOTE

For information on the foreseeable misuse of connected hardware components, please refer to the user manual of the respective device.

1.6 Residual risks

WARNING

Despite compliance with all regulations and the implementation of risk-minimizing measures, not all risks can be completely excluded. Residual risks that exist in connection with the use of the product are listed below.

- Improper operation by untrained personnel can result in harm to the patient.
- Incorrect entry of information in the app, or incorrect assignment of patients or images by the user, may result in a misinterpretation of the latter. The consequences can be an unnecessary treatment or delayed treatment of a skin condition.
- Misuse by the user cannot be ruled out completely despite the provision of information for use.

1.6.1 IT-Security

The following residual risks regarding IT-Security cannot be ruled out completely despite the implementation of risk control measures:

- Accessing and **using another user's credentials, such as username and password (Spoofing)**
- Maliciously changing or modifying persistent data and the alteration of data in transit (Tampering)
- Performing prohibited operations in a system that lacks the ability to trace the operations (Repudiation)
- Reading a file that one was not granted access to, or reading data in transit (Information disclosure)
- Attempting to deny access to valid users, such as by making a web server temporarily unavailable or unusable (Denial of Service)
- Gaining privileged access to resources in order to gain unauthorized access to information or to compromise a system (Elevation of privilege)

Those residual risks may lead to therapeutic patient data being published along with the name of the patient in the worst case.

2 FotoFinder skeen software

2.1 First login

The following steps are carried out in sequence when logging in for the first time:

1. Select the language
2. Set up a WLAN connection
3. First login to the FotoFinder Hub (cf. chapter 2.4 FotoFinder Hub login)
4. Set up a PIN for skeen (cf. chapter 2.3 PIN (Personal Identification Number))

These steps are also required if you logged out of the FotoFinder Hub the last time you used it.

NOTE

Supported Wi-Fi security standards:

- WEP
 - WPA/WPA2 Personal
 - WPA2/WPA3 Personal
 - WPA3 Personal
 - WPA-Enterprise
 - Public Wi-Fi with web browser confirmation (Captive Portal)
-

2.2 Installation und updates

The FotoFinder skeen software is already installed on your device. Updates are carried out automatically via Internet connection (Wi-Fi).

2.3 PIN (Personal Identification Number)

To use your device, you must assign a four-digit PIN when logging in for the first time. The PIN can be freely selected.

You must enter the PIN at the start of each use.

You can change or activate the PIN in the *Settings / System configuration / Security* menu.

2.4 FotoFinder Hub login

Ihr Hub Account

Für die Nutzung von skeen benötigen Sie einen FotoFinder Hub Account!

Haben Sie bereits einen Account?
Loggen Sie sich einfach mit Ihren bestehenden Zugangsdaten in skeen ein.

1 Noch keinen Account? So einfach geht's:

- Besuchen Sie: hub.fotofinder.de
- Erstellen Sie Ihren Account.

Wichtig: Dies ist nur über einen PC oder ein Tablet möglich, NICHT über skeen.

2 HUB Lizenz aktivieren:

- Gehen Sie zu Einstellungen
- Wählen Sie Rechnungsstellung
- Wählen Sie eine passende Lizenz aus oder aktivieren Sie Ihr Lizenzpaket mit einem erworbenen Voucher-Code.

3 Loggen Sie sich mit skeen ein:
Geben Sie Ihre Zugangsdaten ein oder scannen Sie den QR-Code – diesen finden Sie im Hub unter:

- Einstellungen
- Meine Geräte
- Neues Gerät hinzufügen

Jetzt sind Sie startklar!



HUB Instructions
Wir zeigen Ihnen, wie Sie Ihr Gerät in HUB integrieren.
We show you how to integrate your device into HUB.

Your Hub Account

To use skeen, you need a FotoFinder Hub account!

Already have an account?
Simply log in to skeen with your existing credentials.

1 Don't have an account yet? Here's how to create one:

- Visit: hub.fotofinder.de
- Create your account.

Note: This can only be done via a PC or tablet, NOT through skeen.

2 Activate HUB License:

- Go to Settings
- Select Billing
- Choose a suitable license or activate your license package with a purchased voucher code.

3 Log in with skeen:
Enter your login details or scan the QR code - you can find this in the Hub under:

- Settings
- My Devices
- Add New Device

You're all set!



Your photos are stored in the FotoFinder Hub. This cloud solution stores your images and data securely and can be accessed from anywhere. The combination of skeen and Hub automatically synchronises all data, provides access to the integrated AI and serves as an online portal for the analysis and further use of your images. If skeen does not have an online connection, you can also capture pictures in offline mode.

2.5 General operating information

2.5.1.1 Sleep mode

The device switches to sleep mode after a few minutes of inactivity.

There are various ways to reactivate the device:

- press the function button on the side of the device, or
- swipe your finger across the screen from bottom to top, or
- tap the screen twice

2.5.1.2 Software navigation

The software can be easily controlled using gestures. Swipe the screen from left to right with one finger to navigate backwards, for example.

You will also find a back button in every menu.

2.6 Home screen



After connecting to your FotoFinder Hub account, you will see the home screen of the skreen software.

The device and the software are immediately ready for recording.

Fig. 1: Home screen with example preview image

2.7 Menu Bar



You can open the main menu using the menu button at the top left.
The following sub-menus are available:

- My patients
(cf. 2.8)
- Sessions
(cf. 2.11)
- All images
- Settings
(cf. 2.14)
- About FotoFinder
(cf. 2.13)
- FAQ
- Log out

2.8 Patients

Your skreen displays all the patients saved in your FotoFinder Hub account and their images. You can also create new patients with skreen, which are then also synchronised with the Hub.

2.9 Search and select existing patients



Patient search field

1. Tap in the patient search field to open your patient list.
Alternatively, you can also open the patient list via the menu button at the top left in the *My Patients* submenu.





Fig. 2: Example view patient list

2. Select a patient by tapping on the respective line.
You will return to the home screen and the selected patient will be listed at the top.

2.9.1 Create new patient



1. Open the patient list.




2. Tap on *Add New Patient*.

3. Fill in all the required fields in the following window and tap *Add New Patient* again.

The patient has been created and is currently selected.

2.9.2 Remove current patient from selection



1. Open the patient list.

2. Tap on *Unlink* right next to the patient name.

You return to the home screen without a selected patient.

2.10 Creating images

NOTE

The software supports both overview images and micro images. For overview images, remove the magnetic attachment lens from the device. For micro images, the attachment lens with front cap must be attached to the device.

1. Select a patient from the patient list (cf. chapter 2.9 Search and select existing patients).
2. Adjust the required settings. The following buttons are available:

Polarized / Non-polarized



Polarization offers you a special type of light that minimizes reflections on the skin. By pressing this button, you can switch between polarized and non-polarized light. Polarized light is active by default.

Micro image magnification



(15x, 20x or 40x is possible)

If you tap this button, you can choose between the different zoom levels for micro-image capturing.

No zoom is possible for overview images without an attachment lens.

NOTE

-
- A magnification of 20x is required for the analysis with AI score.
 - Micro images with 15x magnification are well suited for trichoscopic examinations.
-

Brightness



When you tap this button, you can choose between three different brightness levels.

3. To take micro images, place the skreen with the attachment lens on the area of skin to be captured and hold it as still as possible.
For overview images, hold the device (without attachment lens) so that you can see the desired image section in the preview window.

4. Press the shutter release button on the handle of the skeen or tap the preview image to release the shutter.

The captured image is displayed.

The following functions are available on the right-hand side of the screen:



AI Score (cf. 2.12)



Save localisation (cf. 2.10.1) or



Delete image

At the bottom of the screen there are two buttons:



5. Tap on Save if you want to capture more micro images. Alternatively, you can use the button on the handle.

The image you have just captured is displayed in small format at the bottom left and you can see the live image again in the preview window.

6. Create further images as described above.
7. Tap on Session overview, or on the preview thumbnail at the bottom left if you want to see an overview of all images created during this session.
8. Tap *Finish* in the *Session overview* if you do not want to add any more images to this session.

This ends the session, and the images are synchronised from skeen to your Hub account. You can also view the images online via your PC in the Hub and use other functions.

2.10.1 Saving the localisation

You can save a localisation for micro images immediately after capturing the image.


Alternatively, you also have this option at a later stage (cf. chapter 2.10.2 Saving the localisation at a later stage).



Immediately after capturing the image, the button shown in the margin appears on the right-hand side of the screen.

1. Tap on it.

This opens the localisation menu:



2. If required, select a different body view by swiping to the right or left.
3. You can use the two-finger zoom to enlarge the view of the localisation mannequin at the required location.
4. Tap on a body part to set it as the localisation. The body part is marked with a red dot and is also listed at the top right.
5. Tap on *Save*. The selected localisation is now saved in the image details during capture.

Fig. 3: Example view of the localisation menu

2.10.2 Saving the localisation at a later stage



1. Open the image for which you want to save a localisation (cf. chapter 2.11 Sessions).
2. Tap the localisation button in the preview window.
3. Continue as described (cf. chapter 2.10.1 Saving the localisation).

2.11 Sessions

You can find an overview of the capture sessions you have already created under *Sessions*.



1. Tap the menu button at the top left.

2. Tap on *Sessions*.

You will see an overview of your previous capture sessions, grouped by day and patient.

3. Tap on a session to open it and you will see all the individual images captured.

You will find the following buttons on the right-hand side of each image captured:



AI Score
(cf. 2.12)



Save localisation
(cf. 2.10.1)



Delete image

2.11.1 Assigning a session to another patient afterwards

If you have assigned an admission session to the wrong patient, you can change this assignment afterwards:

1. Use the menu button at the top left to open the submenu *Sessions*.

You will see an overview of your previous sessions, grouped by day and patient.

2. Tap on the required session.

The session opens.

3. Tap on the pencil icon to the right of the patient selection field.

The patient list opens.

4. Select a patient by tapping on the corresponding line. Alternatively, you can create a new patient (cf. chapter 2.9 Search and select existing patients).

The images from this session are now assigned to this patient.

2.12 The AI Score



The *AI Screening* menu allows you to assess lesions after capturing with Artificial Intelligence. The FotoFinder software uses a Convolutional Neural Network (CNN) algorithm called AI Score. The sensitivity as well as specificity of the algorithm has been proven in a clinical study.

NOTE

Please note that retrieving the AI Score is not available in all countries.

- The AI Score is based on comparisons with images of malignant skin tumors (melanoma, basal cell carcinoma, lentigo maligna, squamous cell carcinoma, actinic keratosis). The Score indicates how similar a lesion is to typical malignant skin tumors.
- The AI Score is not used to assess the malignancy of the examined lesion! It only provides an assessment of whether a lesion is possibly malignant.

NOTE

The AI Score is based on statistics. Therefore, the accuracy of the AI Score cannot be guaranteed and it is intended only as an additional, supportive assessment tool for the physician. The AI Score is not a substitute for the physician's overall clinical diagnosis!

2.12.1 Requesting the AI-Score (AIMEE)

1. Open the relevant micro image. The AI Score is only available for micro images with 20x magnification.



2. Tap on the AI button.

After a short loading process, the AI Score is displayed.




Fig. 4: AI Score Example view

2.12.2 AI Score result information

The AI Score is designed to assess whether a lesion is potentially malignant. This is merely a confidence score of the algorithm, i.e. an assessment of the similarity to malignant lesions. The AI Score is based on comparisons with images of malignant skin tumors (Melanoma, Basal Cell Carcinoma, Lentigo Maligna, Squamous Cell Carcinoma, Actinic Keratosis). The AI Score makes no statement regarding the medical risk and does not assess the malignancy of the examined lesion.

Lesions with a high score should be observed with great attention.

- 0 - 0.49 inconspicuous, follow-up in a reasonable time
 - 0 - 0.2 inconspicuous
 - 0.21 – 0.49 further clarification necessary
- 0.50 - 1.0 conspicuous, should be observed with great attention

2.13 About FotoFinder

In this software section you will find

- the manufacturer's contact details
- the software version
- your device ID
- terms of service
- Terms & Privacy

Explanation of the symbols:

CE Mark




Manufacturer



Country of origin / Date of manufacture




Model number



Displays the Swiss authorized representative:

Johner Medical Schweiz GmbH, Tafelstattstrasse 13a, 6415 Arth, Switzerland



Medical device



Unique Device Identification



Electronic user manual



eIFU indicator

2.14 Settings

You can adjust several functions in the menu *Settings*.

- AI configuration

Here you can choose between

- *Online*: Access to the AI algorithm via the Hub
- *Offline*: Use of the locally installed classification programme (AI algorithm)

- Camera

Here you can change the camera resolution.

- Synchronization

Here you can see when synchronisation with the Hub last took place. You can start a synchronisation at any time (with an active WLAN connection) using the refresh button.

- System configuration

Here you can make settings for the time zone, Wi-Fi, Security (PIN) or automatic system updates, among other things.

FotoFinder	
EU - KONFORMITÄTSERKLÄRUNG EU - DECLARATION OF CONFORMITY	
Hersteller / Manufacturer: Adresse / address:	FotoFinder Systems GmbH Industriestraße 12 84364 Bad Birmbach Deutschland/Germany
Single Registration Number (SRN):	DE-MF-00007084
Benannte Stelle / Notified Body	TÜV SUD Product Service GmbH Ridderstraße 65 80339 München / Munich Germany
Zertifikats-Nr. / Certificate No.	GLO 115802 0002
Wir erklären hiermit in eigener Verantwortung, dass nachstehendes Produkt. <i>We declare under our sole responsibility that the product</i>	
FotoFinder handyscope pro In the following variants FotoFinder handyscope pro, Version: 1.9 FotoFinder skin, Version 1.0	
Zweckbestimmung / Intended Use:	
<p>The FotoFinder handyscope pro app is a mobile application that works in conjunction with the FotoFinder Hub online cloud. The application is designed for patient management, standardized documentation of microscopic images, and to assist in the initial assessment of skin conditions. The FotoFinder handyscope pro enables digital documentation of intact human skin by healthcare professionals. The microscopic images are stored together with the relevant patient data, which makes it possible to visualize changes in lesions during subsequent follow-up examinations of the patient. The FotoFinder application is used in combination with the DermLite handyscope or FotoFinder skin hardware device, which allow to capture microscopic images using a mobile device.</p> <p>The following features are available:</p> <ul style="list-style-type: none">- Acquisition and management of patient data- Capturing and managing microscopic images- Documenting and managing prescriptions- Assigning a session to patient- Assigning a body site (localization) to an image- Requesting a second opinion (Second Opinion) from experts (only for handyscope pro, not for skin)- Request AI score (Artificial Intelligence) <p>The FotoFinder handyscope pro app connects online with the MolesAnalyzer pro algorithms to generate the AI score. These functions are only accessible via paid credits. Credit management is only available through a FotoFinder Hub account. The app data is synchronized, stored and managed via this cloud solution. The FotoFinder handyscope pro is intended for skin lesions. The app must not be used to make or confirm a clinical diagnosis or to monitor the progression of skin cancer. The app may only be used by physicians or healthcare professionals trained in the use of dermoscopy and who have received specific training in the use of the FotoFinder handyscope pro. The AI score is based on statistics. The diagnosis and the prognosis cannot be confirmed by the physician. The application is intended for transient use. In combination with the DermLite handyscope or the FotoFinder skin hardware device, the product is in continuous use for less than 60 minutes during a diagnosis session.</p>	
der Risikoklasse / of risk class:	Ila (Annex VIII MDR)
Basis UDI-DI / Basic UDI-DI:	42601384915A001YV
den Grundlegenden Anforderungen gemäß Anhang I der Medizinprodukteverordnung (EU) 2017/745 entspricht / meets the essential requirements of the regulation (EU) 2017/745.	
Konformitätsbewertungsverfahren / Conformity assessment	(EU) 2017/745, Annex IX Chapters I & III
Diese Erklärung ist gültig, bis sie durch eine neue Version ersetzt wird / This declaration is valid until superseded by a new version.	
Bad Birmbach, 08.10.2024	 Julian Mayer, Authorized Officer